

UNI-DOSE PACKAGING MACHINE



PRODUCT INTRODUCTION



■ UNI-DOSE PACKAGING MACHINE

Automatic subcontracting for oral drugs in inpatient pharmacies and partially agreed prescriptions in outpatient pharmacies. Unicom HIS system obtains the medical order information, automatically packs a single dose of tablets or capsules into the same medicine bag, prints the relevant information of drugs and patients on the medicine bag, and adjusts the quality traceability management.



Rotating Cabinet: Automatically moves the medication cassette for replenishment to the medical staff by the shortest route with light prompts. It also boasts light sensors, halting rotation for safety upon detecting objects.



Intelligent Cassettes: Equipped with 348 intelligent cassettes of three types (S, M, and L), avoiding medication blockage with unique dispensing control technology.



DTA Trays: The two trays, each with 44 compartments, are designed for half and irregularly-shaped drugs, operating alternatively to enhance efficiency. In the case of off-machine drug replenishment, the DTA trays will automatically open, with medication information provided on the screen and replenishment lists printed.



Packaging and Printing: The packaging unit slides out for easy replacement of packaging paper and ink bags. The sealing temperature can be adjusted. Packaging dimensions are 70*70mm, with extended pouch lengths as needed.



12.1-Inch Color Touchscreen: Features a streamlined user interface with Chinese voice prompts, ensuring user-friendly operation.



Built-in Barcode Reader: Supports one-dimensional and two-dimensional barcodes, facilitating the process of medication verification.

Product model:	GP-Q140A	GP-Q160A	GP-Q200A
Drug Storage Type	Drum type	Drum type	Drum type
Intelligent medicine box	120PCS (Big, medium,small)	160PCS (Big, medium,small)	200PCS (Big, medium,small)
DAT drug tray	Standard configuration	Standard configuration	Standard configuration
Dehumidification system	Standard configuration	Standard configuration	Standard configuration
Drug packaging printer	British Domino heat transfer printer	British Domino heat transfer printer	British Domino heat transfer printe
Packer	Standard configuration	Standard configuration	Standard configuration
Smart wristband	Standard configuration	Standard configuration	Standard configuration
Automatic verification and error correction	Standard configuration	Standard configuration	Standard configuration
Touch screen	Standard configuration	Standard configuration	Standard configuration

GP-Q360B GP-Q400B GP-Q450B Product model:

Drug Storage Type	Double layered inner and outer rings type	Double layered inner and outer rings type	Double layered inner and outer rings type
Intelligent medicine box	348PCS (Big, medium,small)	400PCS (Big, medium,small)	440PCS (Big, medium,small)
DAT drug tray	Standard configuration	Standard configuration	Standard configuration
Dehumidification system	Standard configuration	Standard configuration	Standard configuration
Drug packaging printer	British Domino heat transfer printer	British Domino heat transfer printer	British Domino heat transfer printer
Packer	Standard configuration	Standard configuration	Standard configuration
Smart wristband	Standard configuration	Standard configuration	Standard configuration
Automatic verification and error correction	Standard configuration	Standard configuration	Standard configuration
Touch screen	Standard configuration	Standard configuration	Standard configuration





Telephone: LEO: 0086-13391039052; GINA: 0086-18770083375

E-mail : info@g-healthy.com Website : www.g-healthy.com

Address : No.518, Zhongchen Road, Songjiang District, Shanghai, China